Protein A beads

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Table 1 affinity data for Protein A and Protein G for different sources and subtypes of IgG. ++++ = Strong Binding; ++~+++ = Medium Binding; += Weak Binding; +/- = Weak or No Binding; - = No Binding.

Characteristics	Description
Product content	10mg/ml magnetic beads in specific protective buffer
Beads size	~200nm
Magnetization	Superparamagnetic
Coupled protein	Recombinant Protein A
M.W. of protein	~25kDa (Protein A)
Antibody concentration	≥0.6mg Protein A per ml beads
Binding capacity	≥ 0.7mg human IgG per ml beads
Specificity	Antibodies from many different species, including mouse, human, rabbit, cow, goat and sheep
Elution method	Elution with acid, competing peptide or SDS-PAGE loading buffer
Application	IP, Co-IP, Protein purification

Table 2 The main related indicators of Protein A beads.

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Component	K1303-1ml	K1303-5*1ml
Protein A beads	1 ml	5*1 ml

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Problem	Possible Causes	Solution
Very few or no target protein exists in the eluate.	Protein is not completely eluted.	Change elution methods.
	No target protein expressed.	Make sure the protein of interest contains the HA-tag by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate.
		2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Interfering substance is present in sample.	Lysates containing high concentration of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided.
Detection system is inadequate. (WB)	1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity.	
	2. Verify that the transfer was adequate by using prestained protein marker or staining the membrane with Ponceau S.	
	3. Use fresh detection substrate or try a different detection system.	
Background is too high.	Proteins bind nonspecifically to the antibody, insufficient washing on magnetic beads, or the microcentrifuge tubes.	1. Pre-clear lysate with Normal IgG to remove nonspecific binding proteins.
		2. After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before separation.
	Washes are insufficient.	1. Increase the number of washes.
		2. Prolong duration of the washes, incubating each wash for at least 15 minutes.
Multiple protein bands found in the eluate.	The protein is not stable at room temperature.	3. Increase the salt and/or detergent concentrations in the wash solutions.
		4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
	Protein degradation due to proteases activity during purification process.	Purify the target protein at lower temperature, such as 4°C.
Non-specific binding.	Non-specific binding.	Add protease inhibitors to cell lysate.
		1. Prepare cell lysate again.
		2. Add additional wash steps.

Table 3 Other common problems, causes and solutions

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